

LabVIEW Version Compatibility Chart

Scott Hannahs, NHMFL

Saved in LabVIEW Version	Loaded in LabVIEW 2018	Loaded in LabVIEW 2017	Loaded in LabVIEW 2016	Loaded in LabVIEW 2015	Loaded in LabVIEW 2014	Loaded in LabVIEW 2013	Loaded in LabVIEW 2012	Loaded in LabVIEW 2011	Loaded in LabVIEW 2010	Loaded in LabVIEW 2009	Loaded in LabVIEW 8.6	Loaded in LabVIEW 8.5.x	Loaded in LabVIEW 8.2.x	Loaded in LabVIEW 8.0.x	Loaded in LabVIEW 7.1.x	Loaded in LabVIEW 7.0.x	Loaded in LabVIEW 6.1.x	Loaded in LabVIEW 6.0.x	Loaded in LabVIEW 5.1.x	Loaded in LabVIEW 5.0.x	Loaded in LabVIEW 4.0.x	Loaded in LabVIEW 3.0.x	Loaded in LabVIEW 2.5.x	Loaded in LabVIEW 2.2.x	Loaded in LabVIEW 1.2	Loaded in LabVIEW 1.0	
2018	OK	SFP	SFP	SFP	SFP	SFP	SFP	SFP	SFP	SFP	SFP	SFP	SFP	SFP	No	No	No	No	No	No	No	No	No	No	No	No	
2017	OK	OK	SFP	SFP	SFP	SFP	SFP	SFP	SFP	SFP	SFP	SFP	SFP	SFP	No	No	No	No	No	No	No	No	No	No	No	No	
2016	OK	OK	OK	SFP	SFP	SFP	SFP	SFP	SFP	SFP	SFP	SFP	SFP	SFP	No	No	No	No	No	No	No	No	No	No	No	No	
2015	OK	OK	OK	OK	SFP	SFP	SFP	SFP	SFP	SFP	SFP	SFP	SFP	SFP	No	No	No	No	No	No	No	No	No	No	No	No	
2014	OK	OK	OK	OK	OK	SFP	SFP	SFP	SFP	SFP	SFP	SFP	SFP	SFP	No	No	No	No	No	No	No	No	No	No	No	No	
2013	OK	OK	OK	OK	OK	OK	SFP	SFP	SFP	SFP	SFP	SFP	SFP	SFP	No	No	No	No	No	No	No	No	No	No	No	No	
2012	OK	OK	OK	OK	OK	OK	OK	SFP	SFP	SFP	SFP	SFP	SFP	SFP	No	No	No	No	No	No	No	No	No	No	No	No	
2011	OK	OK	OK	OK	OK	OK	OK	OK	SFP	SFP	SFP	SFP	SFP	SFP	No	No	No	No	No	No	No	No	No	No	No	No	
2010	OK	OK	OK	OK	OK	OK	OK	OK	OK	SFP	SFP	SFP	SFP	SFP	No	No	No	No	No	No	No	No	No	No	No	No	
2009	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	SFP	SFP	SFP	SFP	No	No	No	No	No	No	No	No	No	No	No	No	
8.6	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	SFP	SFP	SFP	No	No	No	No	No	No	No	No	No	No	No	
8.5.x	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	SFP	SFP	No	No	No	No	No	No	No	No	No	No	No	
8.2.x	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	SFP	No	No	No	No	No	No	No	No	No	No	No	No	
8.0.x	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	SFP	No	No	No	No	No	No	No	No	No	No	No	
7.1.x	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	SFP	No	No	No	No	No	No	No	No	No	No	
7.0.x	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	SFP	No	No	No	No	No	No	No	No	No	
6.1.x	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	SFP	No	No	No	No	No	No	No	No	
6.0.x	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	SFP	No	No	No	No	No	No	No	
5.1.x	No	No	No	No	No	No	No	No	No	No	No	No	OK	OK	OK	OK	OK	OK	OK	OK	SFP	No	No	No	No	No	
5.0.x	No	No	No	No	No	No	No	No	No	No	No	No	OK	OK	OK	OK	OK	OK	OK	OK	OK	No	No	No	No	No	
4.0.x	No	No	No	No	No	No	No	No	No	No	No	No	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	No	No	No	No	
3.0.x	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	OK	OK	OK	OK	No	No	No	No	
2.5.x	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	OK	OK	No	No	
2.2.x	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	OK (Mac)	No	OK	No	No
1.2.X	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	OK	OK	No	
1.0.X	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	No	OK	OK	

The left hand side gives the version of LabVIEW for the saved VI, the top row gives the version of LabVIEW that can load that VI.

SFP = VI can be opened if "saved for previous" version was used

OK = This version of the VI can be opened by this newer version of LabVIEW

NO = This version of the VI can not be opened or saved for previous version